



Skype: chinacaremed

Facebook: China care medical

Wechat: joseph20130708

Whatsapp: +86 139 2648 0041

Email: joseph@chinacaremedical.com

Website: www.chinacaremedical.com

Add: Room B205, Changsheng Business Building, Daguan South Road No 26, Tianhe, Guangzhou, China

40L, 60L OXYGEN CONCENTRATOR



Model: **CCS-JAY40, CCS-JAY60**

Application:

1. Medical use
2. Industrial use such as welding, glass cutting etc.
3. Fish farming, aquaculture
4. Veterinary use
4. Oxygen supply for ozone generator
5. Diffuse type oxygen supply use in high altitude area or activity center etc.

Specifications

Model: CCS-JAY40 / CCS-JAY60

Flow Rate: 1-40LPM / 1-60LPM

Purity: 93% ($\pm 3\%$)

Outlet Pressure: 0.04-0.07MPa

Sound Level: $\leq 58\text{db}$ / $\leq 60\text{db}$

Power Consumption: $\leq 2200\text{w}$ / $\leq 3300\text{w}$



Skype: chinacaremed

Facebook: China care medical

Wechat: joseph20130708

Whatsapp: +86 139 2648 0041

Email: joseph@chinacaremedical.com

Website: www.chinacaremedical.com

Add: Room B205, Changsheng Business Building, Daguan South Road No 26, Tianhe, Guangzhou, China

LCD Display: Switch times, operating pressure, present working time, accumulated working time

Power Supply: AC220V $\pm$ 10%, 50/60 $\pm$ 1Hz, AC110V $\pm$ 10%, 50/60 $\pm$ 1Hz

Alarm: Low purity alarm, High&low pressure alarm, Power failure alarm

Net Weight: 120Kgs/ 180Kgs

Size: 820 $\times$ 355 $\times$ 1160mm / 820 $\times$ 355 $\times$ 1730mm