



Skype: chinacaremed

Facebook: China care medical

Wechat: joseph20130708

Whatsapp: +86 139 2648 0041

Email: joseph@chinacaremedical.com

Website: www.chinacaremedical.com

Add: Room B205, Changsheng Business Building, Daguan South Road No 26 ,Tianhe,Guangzhou, China

MEDICAL OXYGEN CONCENTRATOR

Model: **CCS-3BW, CCS-5BW**



Optional pulse oximeter



Low maintenance,
easy to care for



Higher outlet pressure
to ensure proper flow
delivery



Skype: chinacaremed

Facebook: China care medical

Wechat: joseph20130708

Whatsapp: +86 139 2648 0041

Email: joseph@chinacaremedical.com

Website: www.chinacaremedical.com

Add: Room B205, Changsheng Business Building, Daguan South Road No 26, Tianhe, Guangzhou, China

Specifications

Model: CCS-3BW/ CCS-5BW

Flow Rate: 0-3LPM / 0-5LPM

Purity: 93% ($\pm 3\%$)

Outlet Pressure: 0.04-0.07MPa

Sound Level: ≤ 43 db

Power Consumption: ≤ 300 w / ≤ 400 w

LCD Display: Switch times, operating pressure, present, working time, accumulated working time, presetting time from 10mins to 40hours

Alarm: Power Failure alarm, High&low pressure alarm

Net Weight: 16Kgs / 23Kgs

Size: 360×300×600mm

Optional

1. Nebulizer (Atomization)
2. Low purity alarm/Low flow alarm
3. High temperature alarm
4. Pulse oximeter
5. Maintenance alarm