



Skype: chinacaremed

Facebook: China care medical

Wechat: joseph20130708

Whatsapp: +86 139 2648 0041

Email: joseph@chinacaremedical.com

Website: www.chinacaremedical.com

Add: Room B205, Changsheng Business Building, Daguan South Road No 26 ,Tianhe,Guangzhou, China

HEMOCARE OXYGEN CONCENTRATOR



Model: CCS-3AW, CCS-5AW



Skype: chinacaremed

Facebook: China care medical

Wechat: joseph20130708

Whatsapp: +86 139 2648 0041

Email: joseph@chinacaremedical.com

Website: www.chinacaremedical.com

Add: Room B205, Changsheng Business Building, Daguan South Road No 26, Tianhe, Guangzhou, China



Low noise,
takes up less space



Innovative big
LCD display



USA imported humidifier,
6PSI pressure release

Specifications

Model: CCS-3AW/ CCS-5AW

Flow Rate: 0-3LPM / 0-5LPM

Purity: 93% ($\pm 3\%$)

Outlet Pressure: 0.04-0.07MPA

Sound Level: $\leq 40\text{db}$

Power Consumption: $\leq 280\text{w}$ / $\leq 300\text{w}$

LCD Display: Switch times, operating pressure, present, working time, accumulated working time, presetting time from 10mins to 40hours

Alarm: Power Failure alarm, High&low pressure alarm

Net Weight: 16Kgs / 16 Kgs

Size: 350×280×540 mm

Optional

1. Nebulizer (Atomization)
2. Low purity alarm/Low flow alarm
3. High temperature alarm
4. Pulse oximeter
5. Maintenance alarm