



**Skype:** chinacaremed

**Facebook:** China care medical

**Wechat:** joseph20130708

**Whatsapp:** +86 139 2648 0041

**Email:** info@chinacaremedical.com

**Website:** www.chinacaremedical.com

**Add:** Room C208, Changsheng Business Building, Daguan South Road No 26, Tianhe, Guangzhou, China

---

## Model CCLTK biological tissue roasting machine



### Performance characteristics:

- 1, the body surface adopts APS electrostatic spraying technology, high strength and corrosion resistance.
- 2, the machine uses single-chip microcomputer to control heating, so that the display temperature is accurate; multiple overheating protection, safe and reliable.
- 3, easy to operate, digital display temperature, one machine multi-purpose.
- 4, the oven is equipped with a plurality of grill grooves, which can be used for oblique insertion or dyeing.
- 5, the optional side cold light source makes the film look clearer and easier.

### Main technical indicators:

1. The heating of the stalls: room temperature to 99 degrees C
2. The heating of the roast slices: room temperature to 99 degrees C
3. The heating of the baking sheet: room temperature to 99 degrees C
- 4, a baking sheet: three dyed frames
5. The power consumption of the whole machine is less than 1KVA
6. Power supply: AC220V + 10%, 50HZ + 1HZ
- 7, the weight of the whole machine: less than 25kg