




Patient Monitor

Model: CCM90N



Feature:

- 12.1"high resolution color TFT display.
- Five standard parameters: ECG, RESP, NIBP, SPO2, 1-TEMP.
- SpO2 Pitch Tone Variation
- Multi-display:standard, large font, trend coaxis, OxyCRG dynamic
- 7-lead ECG waveforms display simultaneous.
- Large tabular and graphic trends information storage
- Built-in rechargeable lithium battery
- Networking capacity and TCP/IP networking platform.
- Optional parameter:
Touch Screen, Recorder, FM,FHR, TOCO,4-IBP, EtCO2

**China Care Medical Equipment
Co., Ltd.**

**Address: Room B101, Zhongbang
Business Building, No. 87 Puxing
Road, Tianhe, Guangzhou,
www.chinacaremedical.com**