



Skype: chinacaremed

Facebook: China care medical

Wechat: joseph20130708

Whatsapp: +86 139 2648 0041

Email: info@chinacaremedical.com

Website: www.chinacaremedical.com

Add: Room C208, Changsheng Business Building, Daguan South Road No 26 ,Tianhe,Guangzhou, China

CCLS1418 Biological tissue dyeing machine



Features:

1. High-definition & sensitive of touch LCD
2. Unique design for dyeing process
3. Flexible transmission, noiseless
4. Intellectualized infall, outfall & drainage system
5. Special dyeing shelf loading with decades of specimen
6. Streaming water cleaning & tissue roasting function

Main Technical Specifications :

- ◆ Reagent containers: 12/16
- ◆ One cleaning bath (the 9th one)
- ◆ One drying bath (the 1st one)



Skype: chinacaremed

Facebook: China care medical

Wechat: joseph20130708

Whatsapp: +86 139 2648 0041

Email: info@chinacaremedical.com

Website: www.chinacaremedical.com

Add: Room C208, Changsheng Business Building, Daguan South Road No 26, Tianhe, Guangzhou, China

- ◆ Processing time of each container: 0s~59 min and 59 s.
- ◆ Capacity of containers: 1 L, specimen process: 72 pc.
- ◆ Working voltage: AC220 V 50Hz ,AC110 V 60Hz
- ◆ Power: 1000 VA
- ◆ Dimensions: 1130 x 460 x 400 mm
- ◆ Weight: 85 kg
- ◆ Unit Dimensions: 1350*450*400 mm,
- ◆ Gross Weight: 110 kg