

Vital Signs Monitor

NC3 is specially designed for spot checking of patient vitals for important parameters such as NIBP, SPO₂, TEMP and Heart Rate. Using infrared Ear temperature measurement technology, the ear thermometer can quickly read the figure of temperature in 2~3 seconds, as well as transfer and display the figure on monitor by the built-in WiFi between ear thermometer and vital signs monitor NC3. It combines with smart and aesthetic design having super li-ion battery and clear screen, totally meets the fastest support for hospital use.

50 groups of data storage,
easy to review

Rapid NIBP measurement technology,
it shorten the measuring time to 20~30
seconds from former 60~120 seconds,
provides PR figure display

Masimo SPO₂ with PI(perfusion index)
figure display

Power on/off :Short time to start up,
it just takes 1 seconds to go into normal
operation interface

6 inch LED back-light LCD display
provides a clear and distinct view

Using contactless pyroelectric infrared
technology to get the ear temperature,
the measurement time is only needed
2 to 3 seconds; effectively prevents the
risk of cross infection.

Provided PR from SpO₂ or NIBP measrement

NIBP start/stop

Retset / Review

Patient Type



High performance



Quick and easy to use, accurate measurement

- Using contactless pyroelectric infrared technology to get the ear temperature, the measurement time is only needed 2 to 3 seconds; effectively prevents the risk of cross infection
- Built-in WIFI module for TEMP function, patient's temperature data is quickly transferred to NC3 by WIFI, as well as the figure is displayed on screen, convenient to measure different patient, it shorten the time and cuts the clinical application cost
- Rapid NIBP measurement technology, it shorten the measuring time to 20~30 seconds from former 60~120 seconds, provides PR figure display
- Quickly switches the relative parameters for different type of patient such as adult, child and neonate
- Optional Masimo SPO2 with PI(perfusion index) figure display, PI provides caregiver with an indication of the reliability of SPO2 measurement

Smart and lightweight, sturdy and durable

- 6 inch LED back-light LCD display provides a clear and distinct view
- The net weight is 1.25kg, small enough and light enough to ward round monitoring
- Various wall mount and trolley movement solutions
- Precision and hermetic design, it is able to withstand 1.2 meters' drop test
- Liquid protection grade gets to IPX1, effectively prevents the damage by liquid

To meet clinical need, excellent vital signs monitor

- 50 groups of data storage, easy to review
- Super li-on battery provides more than 12 hours of cautious monitoring
- Various accessories of NIBP and SPO2 meet the need of different patients including adult, child and neonate



Vital Signs Monitor