



China Care Medical Equipment Co., Limited

Tel: +86 20 62871996

Mob: +86 13711403953

Contact person: Shally Xie

Email: shallychinacare@gmail.com

Web: <http://www.chinacaremedical.com>

Add: DaXin Building 2203C room, DeZheng North Road No.538, YueXiu District, Guangzhou, China

Laptop 4D B/W Ultrasound Machine

Model: CCB-SS10Plus



Feature:

- 4D imaging function
- 15 inch LED screen, pretty, portable and light
- 2 probe connector
- Full digital imaging technology, crystal-clear Image
- PC based, abundant functions
- Powerful image and report management
- Large volume storage image and cineloop
- DICOM display and transfer(optional)
- Compatible with laser/inkjet printers
- Built-in battery(Work 5 hours)

Specification:

- Monitor: 15 inch LED
- Display mode: B, B+B, B+M, M, 4B, 9B
- Image gray scale: 256 level
- Cine loop: Auto 256 frames, Manual, 4B/9B, Auto/Manual M, B/M
- Image storage: >50,000 frames



China Care Medical Equipment Co., Limited

Tel: +86 20 62871996

Mob: +86 13711403953

Contact person: Shally Xie

Email: shallychinacare@gmail.com

Web: <http://www.chinacaremedical.com>

Add: DaXin Building 2203C room, DeZheng North Road No.538, YueXiu District, Guangzhou, China

- Scan angle: 60-160 degrees, depending on probe
- Scan depth: 20mm~240mm
- Dynamic range: ≥ 100 dB, 4-level adjustable
- Image flip: Up/down, left/right, black/ white, 90/180/270°
- Zoom: Continuously, real-time dynamic PIP
- Focus position: Adjustable
- Focal space: Adjustable
- Measurement: Distance, circumference, area, volume, angle, heart rate, time, velocity, residual urine, etc.
- Software packages: General, OB, GYN, cardiology, urinary, small organ
- Notation: Date, time, name, sex, age, doctor, hospital name, full screen notation
- Body mark: 43 types
- Peripheral port: Video, S-video, USB2.0, DICOM3.0, VGA
- Battery volume: 8,000mAh
- Power consumption (MAX): 75VA
- Dimension: 35.5*36.6*7.5 CM
- Net Weight: 5kgs

Standard Configuration:

One Main Unit, One convex probe, One 4D volume probe, One Power adapter.

Optional:

Convex, Linear, Trans-vaginal, Micro-convex, Rectal Probe; DICOM3.0; Trolley; Laser Printer.





China Care Medical Equipment Co., Limited

Tel: +86 20 62871996

Mob: +86 13711403953

Contact person: Shally Xie

Email: shallychinacare@gmail.com

Web: <http://www.chinacaremedical.com>

Add: DaXin Building 2203C room, DeZheng North Road No.538, YueXiu District, Guangzhou, China

