

DispensMate
Bottle-Top Dispenser
CCL-DispensMate

Easy-to-use volume adjustment

Smooth piston design provides effortless dispensing

Adapters allow fitting to standard reagent bottles

Filling tube adjusts easily for different size bottles

Closure cap protects, against reagent contact and is easy to attach/detach, without taking gloves off.



Features

- Excellent chemical resistance, components are made of PTFE, FEP, BSG, PP
- Fully autoclavable at 121
- Four ranges of bottle-top dispenser covering a volume range from 0.5mL to 50mL
- Easy for cleaning and maintenance
- The optional flexible discharge tube with safety handle permits fast and precise dispensing
- Vapor pressure max. 500mbar, viscosity max. 500mm²/s,
- Applied with S40, GL32, GL38, GL25, GL28 sized adapters C, density max. 2.2g/cm³

Accessories

Cat. No.

Description

| | |
|----------|--------------------------------------|
| 17100553 | GL45, Adapter, 45mm |
| 17000069 | S40, Adapter, 45/40mm |
| 17400131 | GL32, Adapter, 45/32mm |
| 17400132 | GL38, Adapter, 45/38mm |
| 17400133 | GL25, Adapter, 32/25mm |
| 17400134 | GL28, Adapter, 32/28mm |
| 17400037 | Reagent Bottle (Brown, 1L, Optional) |
| 17400021 | Filling tube(m) |
| 17100073 | Discharge tube(m) |

Specifications

| Cat. No. | Volume Range mL | Graduation mL | A _± | | CV _≤ | |
|--------------|-----------------|---------------|----------------|-----|-----------------|----|
| | | | % | μL | % | μL |
| 731100017777 | 0.5-5 | 0.1 | 0.5 | 25 | 0.1 | 5 |
| 731100027777 | 1.0-10 | 0.2 | 0.5 | 50 | 0.1 | 10 |
| 731100037777 | 2.5-25 | 0.5 | 0.5 | 125 | 0.1 | 25 |
| 731100047777 | 5.0-50 | 1.0 | 0.5 | 250 | 0.1 | 50 |